


VITAPAN®



Mould Chart

Stand 10-08



VITA



ovoid



tapering

0 13

41,7
9,5
7,9



0 25

42,8
9,9
8,1



0 34

43,3
9,7
8,0



0 35

43,2
9,8
8,0



0 43

43,6
9,4
8,2



0 86

47,1
9,9
8,9



0 97

49,2
11,5
9,0



***0 98**

48,1
11,8
9,0



0 99

52,5
12,5
9,8



***T 36**

43,7
10,1
8,1



T 53

44,6
9,2
8,1



T 56

44,4
10,0
8,3



T 66

44,6
10,2
8,3



T 67

45,1
10,0
8,4



***T 76**

45,9
10,4
8,2



***T 77**

46,2
10,6
8,5



***T 88**

47,1
10,6
8,7



T 98

50,5
10,8
9,3



T 99

48,0
11,1
9,0





rectangular

X 13

41,5
9,3
7,7



X 66

45,9
10,1
8,9



X 77

47,0
10,1
8,5



X 87

46,9
10,1
8,9




X 96

50,0
10,0
9,3



X 99

50,6
12,7
9,4



The systematic division according to mould groups facilitates the selection and classification of the set of teeth according to the physiognomy of the patient. This also facilitates accurate communication between dentist and dental technician.



square

Z 51

44,6
8,8
8,3



Z 61

45,5
8,7
8,2



Z 74

46,7
9,6
8,8



Z 84

46,7
9,4
8,7



Z 85

47,5
9,7
8,6



Z 97

52,2
10,2
9,5




Z 97

52,2

10,2

9,5



All measurements in millimetres

Lower anteriors

Fully anatomic
incisal edges

Abraded
incisal edges

L3

33,4
8,0
5,0



***L5**

35,3
9,2
5,2



L5L

34,2
9,9
5,1



***L7**

36,2
8,7
5,4



L9

36,1
9,0
5,6



***L11**

36,9
9,1
5,7



L13

37,0
8,7
5,6



L15

38,0
12,2
5,7



***L4**

32,7
8,3
4,7



L8

38,9
9,8
5,8



L10

39,0
9,2
5,9



L12

41,0
9,8
6,2






L14

41,7
10,9
6,1



VITAPAN® Combination Chart




Nr.		Nr. 	Nr. 
O 13	41,7	L4 32,7	40C 10S 20E 28,1 29,0 29,4 29,5 30,5 31,0
O 25	42,8	L3 33,4	40C 10S 20E 28,1 29,0 29,4 29,5 30,5 31,0
O 34	43,3	L3 33,4	41C 11S 21E 21L 29,9 30,0 30,5 30,1 31,3 31,9 32,7 30,7
O 35	43,2	L3 33,4	41C 11S 21E 21L 29,9 30,0 30,5 30,1 31,3 31,9 32,7 30,7
O 43	43,6	L5 35,3	41C 11S 22E 22L 29,9 30,0 32,9 31,6 31,3 31,9 33,7 32,1
O 86	47,1	L11 36,9	42C 12S 23E 23L 31,7 31,1 33,5 34,2 33,4 32,7 34,6 34,8
O 97	49,2	L8 38,9	44C 14S 24E 32,8 33,1 34,8 34,6 35,4 36,1
O 98	48,1	L13 37,0	42C 12S 24E 31,7 31,1 34,8 33,4 32,7 36,1
O 99	52,5	L14 41,7	45C 14S 25E 34,8 33,1 36,7 36,3 35,4 37,3
T 36	43,7	L4 32,7	40C 10S 20E 28,1 29,0 29,4 29,5 30,5 31,0
T 53	44,6	L5 35,3	41C 11S 21E 21L 29,9 30,0 30,5 30,1 31,3 31,9 32,7 30,7
T 56	44,4	L5 35,3	41C 11S 22E 22L 29,9 30,0 32,9 31,6 31,3 31,9 33,7 32,1
T 66	44,6	L5 35,3	42C 12S 22E 22L 31,7 31,1 32,9 31,6 33,4 32,7 33,7 32,1
T 67	45,1	L4 32,7	42C 12S 22E 22L 31,7 31,1 32,9 31,6 33,4 32,7 33,7 32,1
T 76	45,9	L7 36,2	42C 12S 23E 23L 31,7 31,1 33,5 34,2 33,4 32,7 34,6 34,8
T 77	46,2	L11 36,9	42C 12S 24E 31,7 31,1 34,8 33,4 32,7 36,1

VITAPAN® Anteriores/Posteriores:

A well-balanced, high-quality assortment with selected tooth moulds to meet the highest aesthetic and economical demands.

The intricate manual layering technique as well as the individual characterization of the VITAPAN® Teeth offer the best guarantee of achieving a natural appearance also in the case of partial and complete dentures. VITAPAN® Anteriores and Posteriores are available in the VITA SYSTEM 3D-MASTER® shades (with the exception of SYNOFORM) and in VITAPAN® classical (without B1). Shade compatibility with all

VITA materials is guaranteed. This makes the VITAPAN® Teeth ideally suitable for the more demanding techniques such as combined and implantborne prosthetics.

Nr.		Nr. 	Nr. 
T 88	47,1	L11 36,9	44C 14S 24E 32,8 33,1 34,8 34,6 35,4 36,1
T 98	50,5	L15 38,0	45C 14S 24E 34,8 33,1 34,8 36,3 35,4 36,1
T 99	48,0	L11 36,9	44C 14S 25E 32,8 33,1 36,7 34,6 35,4 37,3
X 13	41,5	L4 32,7	41C 11S 20E 29,9 30,0 29,4 31,3 31,9 31,0
X 66	45,9	L9 36,1	42C 12S 22E 22L 31,7 31,1 32,9 31,6 33,4 32,7 33,7 32,1
X 77	47,0	L11 36,9	42C 12S 23E 23L 31,7 31,1 33,5 34,2 33,4 32,7 34,6 34,8
X 87	46,9	L11 36,9	44C 14S 24E 32,8 33,1 34,8 34,6 35,4 36,1
X 96	50,0	L10 39,0	44C 14S 25E 32,8 33,1 36,7 34,6 35,4 37,3
X 99	50,6	L15 38,0	45C 14S 25E 34,8 33,1 36,7 36,3 35,4 37,3
Z 51	44,6	L3 33,4	43C 13S 22E 22L 31,0 31,3 32,9 31,6 32,3 33,9 33,7 32,1
Z 61	45,5	L3 33,4	43C 12S 22E 22L 31,0 31,1 32,9 31,6 32,3 32,7 33,7 32,1
Z 74	46,7	L9 36,1	42C 12S 23E 23L 31,7 31,1 33,5 34,2 33,4 32,7 34,6 34,8
Z 84	46,7	L9 36,1	42C 12S 23E 23L 31,7 31,1 33,5 34,2 33,4 32,7 34,6 34,8
Z 85	47,5	L13 37,0	44C 14S 24E 32,8 33,1 34,8 34,6 35,4 36,1
Z 97	52,2	L12 41,0	45C 14S 25E 34,8 33,1 36,7 36,3 35,4 37,3

VITAPAN® Anteriores/Posteriores - a convincing system:

- Translucency, opalescence and fluorescence correspond to that of natural teeth
- The natural surface structure and mould of each individual tooth guarantee the overall harmony of the entire set of teeth
- The VITAPAN®CUSPIFORM and SYNOFORM posteriors are suitable for all recognized setup methods:
CUSPIFORM - with fully anatomical occlusal surfaces, a cuspal inclination of 23-28° and an optimum degree of intercuspation in the centric position SYNOFORM - with abraded occlusal surfaces, a low degree of intercuspation and more narrowly designed moulds
- VITAPAN® Teeth can be individually designed with VITAVMLC.
- Shade stability also in the case of thin wall thicknesses in the veneer technique

Upper and lower posteriors

CUSPIFORM

40 C



****41 C**



****42 C**




43 C



****44 C**



45 C



SYNOFORM

10 S



11 S



12 S



13 S



14 S



SYNOFORM not available in VITA SYSTEM 3D-MASTER® shades.

VITA SYSTEM 3D-MASTER®

With the unique
VITA SYSTEM 3D-MASTER
all natural tooth shades
are systematically
determined and
completely reproduced.


1018E - 1008 (5) SI

VITA shade taking


VITA shade communication

VITA shade reproduction

VITA shade control



VITA

VITA Zahnfabrik is certified, according to the Medical Device Directive.
VITA PHYSIODENS Acrylic Teeth were certified bear the CE mark  0124
and fulfil the requirements of ISO 22112: 2005.
Date of issue of these directions for use: 10/08

VITA Zahnfabrik
H. Rauter GmbH & Co.KG
Postfach 1338
D-79704 Bad Säckingen
Fon +49 (0) 7761/562-0
Fax +49 (0) 7761/562-299
Hotline: +49 (0) 7761/562-222
Fax +49 (0) 7761/562-446
www.vita-zahnfabrik.com
info@vita-zahnfabrik.com